

MANUAL POLARIMETER



The manual Polari meter is ideal for determining concentration Levels of samples for applications in the medical, educational, pharmaceutical, and cosmetic industries etc.

PRODUCT INCLUDED:
100/200mm Cuvettes, Sodium lamp

SPECIFICATIONS:

Model	WXG-4
Measurement Range	-180° ~ +180°
Scale Value	1°
Vernier	0.05°
Magnifier	3x
Light Source	Sodium Lamp
Wavelength	589.44nm
Round Cuvette	100 or 200mm
Power Requirements	220VAC/50Hz
Dimensions	500(L)×135(W)×330(H)mm
Weight	6.5kg