

NDJ Series Rotational Viscometer

NDJ Series Rotational Viscometer is designed for the determination of viscosity and rheological behavior of fluid and semi-fluid. It can also be used to measure the dynamical viscosity and structural viscosity of mineral oil, lubrication oil, adhesives, paints, dyestuff and printing inks, etc.



NDJ-1 Rotational Viscometer



NDJ-4 Rotational Viscometer

NDJ Series Rotational Viscometer		
Model	NDJ-1	NDJ-4
Measurement Range	0.1~1.0*10 ⁵ mPa.s	10~2.0*10 ⁶ mPa.s
Rotor Speed (rpm)	6/12/30/60	0.3/0.6/1.5/3/6/12/30/60
Accuracy	±5% Newtonian Fluid	
Overall Dimensions	445mm*235mm*380mm	
Net Weight	8kg	