

YD-335 Computer Microtome



A newly developed products after YD-202 and YD-1508, it has adopted internationally advanced design. This microtome is current advanced process example, which have CPU inside. So the customer can operate it easily and conveniently.

Essential Features:

- 1)It adopted the advanced import actuation system so that can slice the tissues more precisely ,more reliable and quiet.
- 2)The LCD display show the slice and trimming thickness ,section counter in Chinese and English version.
- 3)It have special function of sample rebounding .it can avoid the damage due to the attrition between sample and the back of knife so that the slice is more smooth and the life of the knife is more longer .
- 4) The hand well can be locked in any position to make the section as safe and convenient as possible
- 5)the waste bath can be disassembled easily
- 6) It has safety alarm system

Technical Data

- 1) Slice thickness range: 0- 100um
- 2) Thickness for repairing the tissue slice: 0 - 100um
- 3) Horizontal displacement of the sample: 28mm
- 4) Vertical dislocation of the sample:60mm
- 5) Sample retraction :12um
- 6) Lowest slice adjusting graduation: 1um
- 7) Tissue slice adjusting graduation: 1um
- 8) Slice precision: 10
- 9) Maximal slice section: 50 x 45mm
- 10)Dimension: 520 x 450 x 300mm
- 11)Weight : 32 kg .

This document was created with Win2PDF available at <http://www.win2pdf.com>.
The unregistered version of Win2PDF is for evaluation or non-commercial use only.
This page will not be added after purchasing Win2PDF.